



- ◆ The Scooter-Lift II Vehicle Lift, Model VSL-900, is ideal for raising and storing a fully-assembled scooter or wheelchair into the sliding side door of minivans, as well as the passenger side door of full-size vans.
- ◆ The unit requires very little arm strength, balance or effort to use. It's so compact that even when you use it in a minivan, there is still enough room for three passengers to sit comfortably in the rear seat. With the scooter removed, or when traveling without the scooter, entry and exit is easy since the door opening remains unobstructed.
- ◆ Bruno's Scooter-Lift II also features an optional installation assembly that eliminates the need to drill holes in the floor of the vehicle. This provides easier removal when it comes time to transfer the lift to another vehicle.
- ◆ The Scooter-Lift II can be used for either three or four-wheel scooters as well as manual-folding or power wheelchairs weighing up to 200 lb/91kg (heavy-duty upgrade available - 225 lb/103 kg). The Scooter-Lift II's hand-held pendant makes this lift easy to use. The Scooter-Lift II is electrically powered from the vehicle's battery and all lift operations are automatic – just help guide the scooter or wheelchair into place.

Optional Bruno Pow'r Pac



- ◆ Easy to use!
- ◆ Convenient!
- ◆ Works in most minivans, SUVs, station wagons & sedans!

Pow'r Pac easily installs in your vehicle – just plug into the vehicle's power outlet!

The complete, independent power supply for your automotive lift!

- ◆ Expert installation and complete, after-sale service is provided by Bruno's worldwide network of dealers.
- ◆ Only Bruno has over 20 different interior/external scooter, wheelchair and powerchair lifts for nearly any sedan, minivan, sport utility vehicle (SUV), full-size van, pickup truck or crossover vehicle (CUV).
- ◆ All Bruno Lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Honda, Hyundai, Saturn, Chrysler, Ford, Toyota, Lexus, Volvo and Volkswagen.

Using the *Bruno Scooter-Lift II* is easy!



Push the **OUT** button to bring the lift above the scooter at the side door; then attach the docking device.



To lift the scooter, push the **UP** button until the wheels clear the van's floor.



Once the scooter is at the proper height, push the **IN** button.



When the scooter is in the van, push the **DOWN** button to set it on the floor.



The Scooter-Lift II also powers a fully-assembled wheelchair into the van effortlessly.



Bruno manufactures a version of the Scooter-Lift II for virtually all minivans with sliding doors, as well as for the passenger side of full-size vans.

Represented By:



Bruno Independent Living Aids, Inc., 1780 Executive Drive,
P.O. Box 84, Oconomowoc, WI 53066
www.bruno.com